


Sample Preparation for 10x Genomics

A How To Guide: Considerations and Best Practices

Bashir Sadet, PhD
Senior Field Application Scientist




Variables important to the success of single cell experiments


- **Objective of the talk:**
 - Identify variables that are crucial for a successful single cell RNA/ATAC-seq experiment
 - +/- feature barcoding
- **Important variables that will be discussed:**
 - Viability
 - Counting
 - Clumps and debris
 - Absence of inhibitors of enzymatic reactions (EDTA, Mg⁺⁺, Ca⁺⁺)
 - Optimization of protocols beforehand

Chromium Single Cell Gene Expression Solution with Feature Barcoding technology

Workflow Overview




- Inputs:
10x Gel Beads, Reagents and single cells in suspension
- Outputs:
Digital gene expression and cell surface protein expression or CRISPR perturbation profiles from every partitioned cell






Single Cell CNV



Single Cell ATAC



Single Cell Gene Expression &
Immune Profiling

Challenging

Most challenging

gDNA

gDNA (in open chromatin)

mRNA

Suboptimal Integrity is a Leading Failure Phenotype

Accurate Cell Recovery

Goal: +/- 20% of target

Challenging sample: <80% of target

Root Causes:

*Sample integrity

*Technical cell counting

High library complexity (signal)

Goal: Cell type dependent

Challenging sample: Very low value

Root Causes:

*Sample integrity

*RNase

Low background (noise)

Goal >70%, Ideal >90%

Challenging sample: <70%

Root causes:

*Sample integrity

*Insufficient wash

Cell Ranger · neuron_10k_v3 · Cells from a combined cortex, hippocampus and subventricular zone of an E18 mouse

SUMMARY ANALYSIS

Estimated Number of Cells

11,843

Mean Reads per Cell

30,153


Median Genes per Cell

2,679

Sequencing

Number of Reads	357,111,595
Valid Barcodes	97.1%
Sequencing Saturation	49.2%
Q30 Bases in Barcode	97.1%
Q30 Bases in RNA Read	95.4%
Q30 Bases in Sample Index	93.8%
Q30 Bases in UMI	97.0%


Cells



Estimated Number of Cells	11,843
Fraction Reads in Cells	90.8%
Mean Reads per Cell	30,153
Median Genes per Cell	2,679
Total Genes Detected	21,080
Median UMI Counts per Cell	6,570


Quality is Critical

Clean



- Aggregates/clumps
- Subcellular debris
- Free-floating RNA/DNA

Healthy



- Biological decomposition
- RNA leakage (background)
- RNA degradation (signal)

Intact



- Physical decomposition
- RNA leakage (background)
- RNA degradation (signal)

Goal is to minimize

General Cell Handling Recommendations

Debris/Aggregate Removal

- Use a cell strainer to remove aggregates or debris from washed cells
- The presence of cell aggregates, debris and/or fibers can result in inaccurate cell counts
- GEM generation occurs in microfluidic channels that are narrower than the typical human hair (i.e. $< 100\ \mu\text{m}$) and the presence of cell debris or large aggregates may clog or wet the chip



Cell Counting

- Quantitate cells accurately before loading into the system
 - Approximately 65% loaded cells will be recovered
 - To maximize the likelihood of achieving the desired recovery target, the optimal input cell concentration is 700-1200 cells/ μl
 - Recommended range: 500 to 10,000 recovered cells
 - Under- or over-loading may impact application performance



Storage of Single Cell Suspensions

- Cell suspensions should always be kept on ice and where possible proceed with cell loading immediately after sample preparation
 - Ideally incubation time should be kept to a minimum ($< 30\ \text{min}$)
- Some cell types are more fragile and cell viability may decrease significantly if not processed and loaded immediately




Cell Washing

Washing isolated cells


- Transfer cells in media to a 2 mL Eppendorf tube
- Spin down cells to form pellet
 - Depending on cell size and concentration, pellet size varies
- Remove supernatant
- Gently add 1x PBS + 0.04% BSA away from cell pellet
- Gently pipette mix with Wide Bore pipette tip
- Repeat the wash one more time
- Spin down cells to form pellet
- Remove supernatant
- Resuspend cells in 1x PBS + 0.04% BSA with gentle pipette mix
 - **For accurate cell counting, do not invert tubes**
- Adjust to desired cell concentration

Note: PBS can be replaced with most common cell culture buffers and media if cells are unstable in PBS

Jurkat: PBS Washes



Live/Dead Staining of Final Suspension





Spin Speed Optimized to Recover Live Cells


Clear Cell Renal Carcinoma

Colorectal Cancer

Before



After filtering and
dead cell removal
spins




Alternative Buffer and Media

Tested in-house

- Tested input volume: 2.5 and 33 μ l
- **Alternative Buffer:** no influence on performance
 - Dulbecco's Phosphate-Buffered Saline (DPBS)
 - Hank's Balanced Salt Solution (HBSS)
- **Alternative Media:** minimal reduction to no loss in performance
 - Eagle's Minimum Essential Medium (EMEM) + 10% FBS
 - Dulbecco's Modified Eagle Medium (DMEM) + 10% FBS
 - Iscove's Modified Eagle Medium (IMEM) + 10% FBS
 - Roswell Park Memorial Institute (RPMI) + 10% FBS
 - Ham's F12 + 10% FBS
 - 1:1 DMEM/F12 +10% FBS
 - M199




Sample Prep Workflow: Rare Cell Population



Example: Frozen Tissue GEX

Homogenization by Pestle



Protocol goal –

Maintain integrity by working quickly:

- Mechanical dissociation
- Use FACS to remove debris and ambient mRNA
- Collect nuclei in GEM-RT mix and process immediately

Why is this sample challenging?


- Small piece of tissue
- Frozen tissue

References

- 10x Genomics Customer Developed Protocol from Luciano Marletotto

Chromium Single Cell Gene Expression Solution with Feature Barcoding technology

Single-use microfluidics chip



- Up to 8 channels processed in parallel
- 500 to 10000 cells per channel
- 18 minute run time (per chip)
- Up to 30 μm cell diameter tested
- Up to 65 % cell processing efficiency
- Temperature Range 18-28°C

User controlled trade-off between cell numbers and doublet rate


Number of Cells Loaded	Number of Cells	Expected Multiplet Rate (%)*
~800	~500	~0.4
~1,600	~1,000	~0.8
~8,000	~5,000	~3.9
~16,000	~10,000	~7.6

Firmware 4.00 or higher is required for running this assay

Single Cell Gene Expression Workflow

Determine Targeted Number of Cells

- Cell Suspension Volume Calculator Table is provided in the User Guide
- Targeted Cell Recovery numbers take into account a 65% cell processing efficiency
- For Example: Targeting 5000 cells with a cell concentration of 1000 cells/ μ l:
 - Add 35.0 μ l water to aliquoted Master Mix
 - Add 8.3 μ l cell suspension to Master Mix containing water
- Remember to resuspend the cells gently when taking aliquots to count or add to Master Mix



Cell Suspension Volume Calculator Table

(for step 1.2 of Chromium Next GEM Single Cell 3' v3.1 protocol)

Volume of Cell Suspension Stock per reaction (μ l) | Volume of Nuclease-free Water per reaction (μ l)

Cell Stock Concentration (Cells/ μ l)	Targeted Cell Recovery										
	500	1000	2000	3000	4000	5000	6000	7000	8000	9000	10000
100	8.3 35.0	16.5 26.7	33.0 10.2	n/a	n/a	n/a	n/a	n/a	n/a	n/a	n/a
200	4.1 39.1	8.3 35.0	16.5 26.7	24.8 18.5	33.0 10.2	41.3 2.0	n/a	n/a	n/a	n/a	n/a
300	2.8 40.5	5.5 37.7	11.0 32.2	16.5 26.7	22.0 21.2	27.5 15.7	33.0 10.2	38.5 4.7	44.0 -0.8	n/a	n/a
400	2.1 41.1	4.1 39.1	8.3 35.0	12.4 30.8	16.5 26.7	20.6 21.6	24.8 18.5	28.9 14.3	33.0 10.2	37.1 6.1	41.3 2.0
500	1.7 41.6	3.3 39.9	6.6 36.6	9.9 33.3	13.2 30.0	16.5 26.7	19.8 23.4	23.1 20.1	26.4 16.8	29.7 13.5	33.0 10.2
600	1.4 41.8	2.8 40.5	5.5 37.7	8.3 35.0	11.0 32.2	13.8 27.5	16.5 26.7	19.3 24.0	22.0 21.2	24.8 18.5	27.5 15.7
700	1.2 42.0	2.4 40.8	4.7 38.5	7.1 36.1	9.4 33.8	11.8 31.4	14.1 29.1	16.5 26.7	18.9 24.3	21.2 22.0	23.6 19.6
800	1.0 42.2	2.1 41.1	4.1 39.1	6.2 37.0	8.3 35.0	10.3 31.9	12.4 30.8	14.4 28.8	16.5 26.7	18.6 24.6	20.6 22.6
900	0.9 42.3	1.8 41.4	3.7 39.5	5.5 37.7	7.3 35.9	9.2 34.0	11.0 32.2	12.8 30.4	14.7 28.5	16.5 26.7	18.3 24.9
1000	0.8 42.4	1.7 41.6	3.3 39.9	5.0 38.3	6.6 36.6	8.3 35.0	9.9 33.3	11.6 31.7	13.2 30.0	14.9 28.4	16.5 26.7
1100	0.8 42.5	1.5 41.7	3.0 40.2	4.5 38.7	6.0 37.2	7.5 35.7	9.0 34.2	10.5 32.7	12.0 31.2	13.5 29.7	15.0 28.2
1200	0.7 42.5	1.4 41.8	2.8 40.5	4.1 39.1	5.5 37.7	6.9 36.3	8.3 35.0	9.6 33.6	11.0 32.2	12.4 30.8	13.8 29.5
1300	0.6 42.6	1.3 41.9	2.5 40.7	3.8 39.4	5.1 38.1	6.3 36.9	7.6 35.6	8.9 34.3	10.2 33.0	11.4 31.8	12.7 30.5
1400	0.6 42.6	1.2 42.0	2.4 40.8	3.5 39.7	4.7 38.5	5.9 37.3	7.1 36.1	8.3 35.0	9.4 33.8	10.6 32.6	11.8 31.4
1500	0.6 42.7	1.1 42.1	2.2 41.0	3.3 39.9	4.4 38.8	5.5 37.7	6.6 36.6	7.7 35.5	8.8 34.4	9.9 33.3	11.0 32.2
1600	0.5 42.7	1.0 42.2	2.1 41.1	3.1 40.1	4.1 39.1	5.2 38.0	6.2 37.0	7.2 36.0	8.3 35.0	9.3 33.9	10.3 32.9
1700	0.5 42.7	1.0 42.2	1.9 41.3	2.9 40.3	3.9 39.3	4.9 38.3	5.8 37.4	6.8 36.4	7.8 35.4	8.7 34.5	9.7 33.5
1800	0.5 42.7	0.9 42.3	1.8 41.4	2.8 40.5	3.7 39.5	4.6 38.6	5.5 37.7	6.4 36.8	7.3 35.9	8.3 35.0	9.2 34.0
1900	0.4 42.8	0.9 42.3	1.7 41.5	2.6 40.6	3.5 39.7	4.3 38.9	5.2 38.0	6.1 37.1	6.9 36.3	7.8 35.4	8.7 34.5
2000	0.4 42.8	0.8 42.4	1.7 41.6	2.5 40.7	3.3 39.9	4.1 39.1	5.0 38.3	5.8 37.4	6.6 36.6	7.4 35.8	8.3 35.0


Single Cell Gene Expression Workflow

cDNA Amplification

- Prepare cDNA Amplification Reaction Mix
 - Primers must be selected based on whether generating only Gene Expression Libraries, or Gene Expression and Cell Surface Protein or CRISPR Libraries.
- Add 65 µl cDNA Amplification Reaction Mix to 35 µl of sample

Lid Temperature	Reaction Volume	Run Time
105°C	100 µl	~30-45 min
Step	Temperature	Time
1	98°C	00:03:00
2	98°C	00:00:15
3	63°C Version Specific Updated Temperature	00:00:20
4	72°C	00:01:00
5	Go to Step 2, see table below for total # of cycles	
6	72°C	00:01:00
7	4°C	Hold


- Use Cell Load to determine number of Total Cycles for cDNA amplification



Cell Load	Total Cycles
<500	13
500–6,000	12
>6,000	11

Sequencing Depth for Typical Samples

Gene Expression Libraries




- 20,000 raw read pairs per cell is the recommended minimum sequencing depth for typical samples.
- Given variability in cell counting and loading, extra sequencing may be required if the cell count is higher than anticipated.

Single Cell Gene Expression Workflow

Carrying forward only 25% of Total cDNA is Sufficient to Capture Full Library Complexity

Correlation: cDNA Input vs Library Complexity




Library Construction Input Mass and SI PCR Cycle Examples


Cell Type	Targeted Cell Recovery	Total cDNA Yield (ng)	cDNA Input into Fragmentation		SI PCR Cycle Number
			Volume (µl)	Mass (ng)	
High RNA Content	Low	250 ng	10 µl	62.5 ng	13
	High	1900 ng	10 µl	475 ng	10
Low RNA Content	Low	1 ng	10 µl	0.25 ng	16
	High	200 ng	10 µl	50 ng	12

Single Cell Gene Expression Workflow

Visually Inspect Gems to Confirm the Presence of a Uniform Emulsion



All liquid levels are similar in volume and opacity without air trapped in the pipette tips.




Pipette tips C and E indicate a wetting failure. Pipette tip C contains partially emulsified GEMs. Emulsion is absent in pipette tip E. Pipette tip G indicates a reagent clog.

QC Metrics: Cell Calling


Wetting Failure

- A comparison of cell calling between a typical sample and a sample with wetting failure.
- Estimated number of cells is unreliable in such failure modes.

Typical Sample



Loss of Single Cell Behavior (eg. Wetting failure)




Algorithm has trouble discerning cells from the background

QC Metrics: Cell Calling


Clog

- A comparison of cell calling between a typical sample and a sample with a clog.
- Similar barcode rank profile can be seen for example for samples with clogs, low sequencing depth.
- The number of cells recovered will likely be lower than the targeted cell numbers. But the data for the cells called is usable.


Typical Sample



Low Barcode counts
(eg. Clog, low-depth)



Reveal Antigen Specificity with Feature Barcoding Technology



Antibody Preparation

Custom Conjugation

Specific Reagents and Consumables

Vendor	Item	Part Number
Expedeon	ThunderLink Plus	425-0300
IDT	Feature Barcode Oligonucleotide (see Table 1 & Appendix)	-
	100 µg Purified Azide-free Antibody	

To perform custom conjugation, a customer needs to acquire:

- ThunderLink Plus kit
- Feature Barcode Oligonucleotide
 - 10 nmoles, HPLC purified, lyophilized
 - NOT** resuspended in TE buffer
 - The BC structure must be compatible with desired assay
 - The BC sequence should be selected from the provided whitelist
 - This ensures compatibility if pooled with commercially available Feature Barcoding technology compatible antibodies
- 100 µg Antibody

10x Genomics Protocol	Feature Barcode Oligonucleotide Sequence			
Single Cell 3' v3 – Cell Surface Protein (CG000185)	/5AmMC12/GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCTNNNNNNNNNN-NNNNNNNNNNNNNNNN-NNNNNNNNNGCTTTAAGGCCGGTCCTAGCAA	TruSeq Read 2	10 nt	Feature Barcode (15 nt) 9 nt Capture Sequence 1
Single Cell V(D)J – Cell Surface Protein (CG000186)	/5AmMC12/CGGAGATGTGTATAAGAGACAGNNNNNNNNNN-NNNNNNNNNNNNNNNN-NNNNNNNNCCCATATAAGAAA	Nextera partial Read 2	10 nt	Feature Barcode (15 nt) 9 nt Capture Sequence

Note: Consult Barcode Whitelist for Custom Feature Barcoding conjugates (Document CG000193), for more information

Dissociation Example: Single Cell ATAC-seq

Fresh mouse breast

Mince with scalpel
blade on ice

Lysis buffer
(1x detergent)

Lysis buffer
(0.1x detergent)



10 min on ice

Spin 5 min 500 rcf


Wash

Resuspend in diluted
nuclei buffer

Signal-to-noise




Sensitivity



Results: **Gentler lysis**
produced higher
quality ATAC data

Sample Quality Assessment

Nuclei Quality



60x Magnification/Brightfield


A: High-quality nuclei have well-resolved edges. Optimal quality for single cell ATAC libraries.

B: Mostly intact nuclei with minor evidence of blebbing. Quality single cell ATAC libraries can still be produced.


C: Nuclei with strong evidence of blebbing. ***Proceed at your own risk.***

D: Nuclei are no longer intact. ***Do not proceed!***


Filter Barcodes



- Some barcodes have very few *overlaps with called peaks*. Should they be kept?



Ideal library



- High number of neutrophils and dead cells
- Low in targeting/ ATAC peak signal/ TSS score
- Insert size distribution is enriched in mononucleosome fragments, no DNA twist

Summary of Key Lessons

When our standard guidance isn't applicable:

- Treat cells gently and minimize decomposition
 - Use gentl(er) lysis conditions
 - Reduce wash steps*
 - Use a swinging bucket centrifuge
 - Keep cells in media + FBS instead of PBS**
- Work quickly
 - Consider sorting, it is a versatile tool for sample prep
 - Minimize unnecessary handling steps
- *Consider the benefits and drawbacks of every different technique*

*Cell surface protein analysis requires thorough washing

**ATAC has specific buffer formulation

Sample Prep Support




10x-pert Webinars

High level overviews of Sample Prep




Demonstrated Protocols & Technical Notes

Learn foundational techniques



Customer Developed Protocols & Publications

Find specific examples



World-Class Support Team

Talk to the experts



Resources

Demonstrated Protocols




User Guides



Application Notes




Public Data Sets



How-To Videos



Scientific Seminars



<https://support.10xgenomics.com/>
<https://www.10xgenomics.com/10x-university/>

10x Support Contact Information

Support Overview:

<http://support.10xgenomics.com>

Q&A Knowledgebase:

<https://kb.10xgenomics.com/hc/en-us>

Please send questions, comments, and feedback to:

support@10xgenomics.com